

Target Materials

	SC7620 Mini Sputter Coater	SC7640 SEM Sputter Coater Auto/Manual	K500X & K550X Sputter Coaters	K650X & K650XT Large Sputter Coater 6" (3 targets required)	SC7680 Large Chamber Sputter Coater	K575X Turbo Chromium Coater	K575XD Turbo Dual Head Coater	K675X Turbo Large Chromium Coater 8" (3 targets required)	K675XD Dual Head Turbo Chromium Coater	SC3000 Turbo Large Chromium Coater 8" (3 targets required)
Aluminium (1.0mm thick)	N/A	N/A	N/A	N/A	N/A	TK8875	TK8875	TK8875	TK8875	TK8875
Carbon (1.5mm thick)	N/A	N/A	N/A	N/A	N/A	TK8869	TK8869	TK8869	TK8869	TK8869
Chromium (0.3mm thick)**	N/A	N/A	N/A	N/A	N/A	TK8845	TK8845	TK8845	TK8845	TK8845
Cobalt (0.1mm thick)	N/A	N/A	N/A	N/A	N/A	TK8900	TK8900	TK8900	TK8900	TK8900
Copper (0.1mm thick)	SC502-314H	SC510-314CU	N/A	N/A	SC510-314CU	SC502-314H	SC502-314H	SC502-314H	SC502-314H	SC502-314H
Gold (0.1mm thick)**	SC502-314A	SC510-314A	TK8842	TK8842	SC510-314A	SC502-314A	SC502-314A	SC502-314A	SC502-314A	SC502-314A
Gold/Palladium (0.1mm thick)	SC502-314B	SC510-314B	TK8843	TK8843	SC510-314B	SC502-314B	SC502-314B	SC502-314B	SC502-314B	SC502-314B
Indium Tin Oxide (9mm thick)	N/A	N/A	N/A	N/A	N/A	TK8907	TK8907	TK8907	TK8907	TK8907
Iridium (0.3mm thick)	N/A	N/A	N/A	N/A	N/A	TK8899	TK8899	TK8899	TK8899	TK8899
Iron (0.1mm thick)	N/A	N/A	N/A	N/A	N/A	TK8897	TK8897	TK8897	TK8897	TK8897
Magnesium (0.3mm thick)	N/A	N/A	N/A	N/A	N/A	TK8905	TK8905	TK8905	TK8905	TK8905
Molybdenum (0.1mm thick)	N/A	N/A	N/A	N/A	N/A	TK8903	TK8903	TK8903	TK8903	TK8903
Nickel (0.1mm thick)	N/A	SC510-314D	TK8874	TK8874	SC510-314D	SC502-314D	SC502-314D	SC502-314D	SC502-314D	SC502-314D
Palladium (0.1mm thick)	SC502-314G	SC510-314G	TK8885	TK8885	SC510-314G	SC502-314G	SC502-314G	SC502-314G	SC502-314G	SC502-314G
Platinum (0.1mm thick)**	SC502-314C	SC510-314E	TK8844	TK8844	SC510-314C	SC502-314C	SC502-314C	SC502-314C	SC502-314C	SC502-314C
Platinum/Palladium (0.1mm thick)	TK8898	N/A	TK8898	TK8898	N/A	TK8878	TK8878	TK8878	TK8878	TK8878
Silver (0.1mm thick)	SC502-314E	SC510-314E	TK8873	TK8873	N/A	SC502-314E	SC502-314E	SC502-314E	SC502-314E	SC502-314E
Tantalum (0.1mm thick)	N/A	N/A	N/A	N/A	N/A	TK8906	TK8906	TK8906	TK8906	TK8906
Titanium (1.5mm thick)**	N/A	N/A	N/A	N/A	N/A	TK8879	TK8879	TK8879	TK8879	TK8879
Tin (0.1mm thick)	N/A	N/A	N/A	N/A	N/A	TK8902	TK8902	TK8902	TK8902	TK8902
Tungsten (0.5mm thick)**	N/A	N/A	N/A	N/A	N/A	TK8846	TK8846	TK8846	TK8846	TK8846

** Alternative thicknesses are available. Please check for availability.

Carbon Supplies

A0830A	Carbon rods 6.15mm dia. x 100mm length	Pack 10
A0832A	Carbon rods 6.15mm dia. x 50mm length (shaped)	Pack 10
A0834A	Carbon rods 3.05mm dia. x 50mm length (shaped)	Pack 10
C5422	Carbon rods 3.05mm dia. x 300mm length	Pack 10
A0819	Carbon fibre (woven cord) - high purity	1m
A0819/5	Carbon fibre (woven cord) - high Purity	5m in 1m lengths
C5421	Carbon fibre (woven cord) - standard purity	1m
C5421-10	Carbon fibre (woven cord) - standard purity	10m
S8650	Hand-held carbon rod shaper (6mm only)	Each
S8651	Hand-held carbon rod shaper (3mm)	Each
SC7605	Motorised carbon rod grinder (state 3mm or 6mm)	Each

