



E5175 SPUTTER TARGETS

Quick Overview

These targets are fitted to the E5175 Sputter Coater. Targets are annular (115mm OD, 92mm ID), 0.1mm thick and supplied as an assembly mounted on an aluminium ring holder.



| | | |
|-----------------------------|------------|------------------------|
| Ordering Information | E5175-314A | Gold (Au) |
| | E5175-314B | Gold/palladium (Au/Pd) |
| | E5175-314C | Platinum (Pt) |
| | E5175-314D | Nickel (Ni) |
| | E5175-314E | Silver (Ag) |

Quorum Technologies Ltd
South Stour Avenue
Ashford
Kent
TN23 7RS
UK
Email:
Website:

Tel: +44 (0) 1233 646332
Fax: +44 (0) 1233 640744
sales@quorumtech.com
http://www.quorumtech.com

For further information regarding any of the other products designed and manufactured by Quorum Technologies, contact your local representative or directly to Quorum Technologies at the address above.