



## SC650 SPUTTER TARGETS

### Quick Overview

These targets are fitted to the SC650 and SC650A sputter coaters. Targets are 0.1mm thick and mounted on aluminium backing plates (set of four).



<b>Ordering Information</b>	SC650-314A	Gold (Au)
	SC650-314B	Gold/palladium (Au/Pd)
	SC650-314C	Platinum (Pt)
	SC650-314D	Nickel (Ni)
	SC650-314E	Silver (Ag)

Quorum Technologies Ltd  
South Stour Avenue  
Ashford  
Kent  
TN23 7RS  
UK  
Email:  
Website:

Tel: +44 (0) 1233 646332  
Fax: +44 (0) 1233 640744  
sales@quorumtech.com  
<http://www.quorumtech.com>

For further information regarding any of the other products designed and manufactured by Quorum Technologies, contact your local representative or directly to Quorum Technologies at the address above.