

EU Declaration of Conformity

Manufacturer:

Quorum Technologies Ltd. Registered Office: Unit 19, Charlwoods Road, East Grinstead, West Sussex, RH19 2HL, United Kingdom. Telephone enquiries to: +44(0)1273 815340. Email: sales@quorumtech.com. Company No: 04273003

The Directives covered by this Declaration:

Simple Pressure Vessels Directive 87/404/EEC
European Pressure Equipment Directive (PED) 97/23/EC (29 November 1999 - SI 1999/2001)

The Products Covered by this Declaration:

E3000/3100 Critical Point Drier

The basis on which conformity is being declared:

The manufacturer hereby declares under his sole responsibility that the products identified above have been tested for compliance with the above directives. The pressure vessel of the products includes no welded element and accordingly the product falls outside the scope of the Simple Pressure Vessels Directive 87/404/EEC. The products were examined against the requirements of European Pressure Equipment Directive (PED) 97/23/EC (29 November 1999). It was determined that the products are outside the scope of the categorisation system as a result of the pressure/volume ratio and therefore are considered to be covered by Sound Engineering Practise. This SEP includes:

- pressure testing of every new product to 2000 psi / 138 bar (working pressure is 1200 psi / 82.7bar)
- the inclusion of an over pressure bursting disc safety device
- submission of sample product to recognised third party organisation for pressure test

The technical documentation required to demonstrate that the products have been examined and comply with SEP are available for inspection by any relevant authority. Categorisation as SEP specifically prohibits application of a CE mark.

Signed: 

Authority: ..Managing Director

Date: 5th November 2001

Attention!

The attention of the specifier, purchaser, installer, or user is drawn to special measures and limitations to use which must be observed when these products are taken into service to maintain compliance with the above directives. Details of these special measures and limitations to use are available on request and are also contained in the product manuals.